


Particle Analysis at the Touch of a Button



Litesizer™ 500



Particle systems can be complex ...

The size and stability of nanoparticles and microparticles are crucial to their performance and function, as well as to their processing and transport properties.

Anton Paar, a leading developer and manufacturer of high performance analytical instruments, has combined its physics and engineering expertise with modern software creativity to create an intuitive particle analyzer that is a joy to use:

The Litesizer™ 500 measures particle size, zeta potential and molecular mass with the dynamic light-scattering (DLS) technology using ingeniously simple software.

It enables rapid and accurate insight into your particle systems and provides the tools for optimizing them by revealing how they change with time, pH, temperature and concentration.



Make
the light
choice

... but measuring them doesn't have to be

The **Litesizer™ 500** will free up your time to concentrate on what your particles are doing, rather than trying to figure out how to use the instrument.




See everything on one page

All input parameters, measurements and analyses are on one page.

Keep track of your samples

The Litesizer™ 500 provides audit trails and compliance with pharmaceutical regulations. Customized reports can be generated in a few seconds and passed on for signing, either electronically or by hand.



See what your particles are doing

A series of measurements lets you see how your particles change with time, temperature, pH or concentration. Results are clearly displayed in different colors so that trends can be identified, while all important values and parameters are logically tabulated below.

Save time

Short start-up times, a one-page workflow, intelligent measurement series, automatic measurement modes, and the fastest zeta potential measurements:
The Litesizer™ 500 saves you time.



Litesizer 500

Technical highlights

Optical bench


The optical bench is the heart of the Litesizer™ 500. Highly sensitive measurement optics enable the accurate detection of even low-intensity signals, while the robust casing reduces the effects of vibrations and ensures that measurements remain unaffected by dust or temperature fluctuations.

One instrument – three detection angles

Choose from back-, side- and forward scattering, or let the Litesizer™ 500 select the best angle for your sample.

Patented ELS technology: cmPALS

The Litesizer™ 500 uses cmPALS, a novel patented PALS technology (European Patent 2 735 870) that defines a new state of the art in ELS optics. The result: zeta potential measurements with the highest accuracy and shortest measurement time.



Continuous transmittance measurements

Continuously measuring light transmitting through the sample enables the device to automatically adjust parameters like focus position, measurement angle, and measurement duration.

Unprecedented size resolution with DLS

The Litesizer™ 500 can precisely resolve bimodal and even trimodal particle mixtures.

The Omega cuvette

The zeta-potential cell has an inverted Ω -shaped capillary tube. This facilitates a homogeneous electric field within the measuring channel, guaranteeing stable and reproducible results.

General specifications:

Temp. control range	0 °C to 90 °C
Light source	Semiconductor laser / 40 mW, 658 nm
Operating temp.	10 °C to 35 °C
Humidity	35 % to 80 % non-condensing
Dimensions (WxDxH)	460 mm x 485 mm x 135 mm
Weight	approx. 18 kg (40 lbs)

Measurement principles

Particle size measured by DLS

Particles suspended in a liquid are constantly undergoing random motion, and the size of the particles directly affects their speed, with smaller particles moving faster than larger ones. In DLS, light passes through the sample, and the scattered light is detected and recorded at a certain angle. The time dependence of the scattered intensity reveals how fast the particles are moving. From this information, it is possible to calculate the average size of the particles as well as their size distribution.

Your benefit:

The Litesizer™ 500 enables accurate and precise particle size measurements as well as easily analyzing the effect of time, pH, temperature and concentration on the particle size. The Litesizer™ 500 provides advanced measurement algorithms which enable you to resolve several different particle sizes in a single suspension.

Particle size specifications

Measuring range	0.3 nm to 10 μm^*
Sensitivity	0.1 mg/mL (lysozyme)
Max. sample concentration	40 % w/v
Accuracy	Better than ± 2 % on certified reference material
Repeatability	Better than ± 2 % on certified reference material
Min. sample volume	20 μL
Measurement angles	15°, 90°, 175°

Molecular mass measured by SLS

The intensity of the scattered light is directly related to the molecular mass of the sample. If the scattering intensity is measured at several different concentrations, then a Debye plot can be generated, the intercept of which provides the molecular weight.


Your benefit:



SLS measurements are simple, fast and non-invasive. They also give you the second virial coefficient, which reflects protein solubility.

Molecular mass specifications

Measuring range	1 kDa to 20 MDa
Sensitivity	0.1 mg/mL (lysozyme)
Repeatability	± 5 %
Measurement angle	90°

* under laboratory conditions





Zeta potential measured by ELS

In electrophoretic light scattering (ELS) the speed of the particles is measured in the presence of an electric field. The faster the particles move, the higher the zeta potential of the particles. In general, a greater-magnitude zeta potential means that the particles will repel each other more strongly, giving a more stable suspension.

Your benefit:

The Litesizer™ 500 uses a patented technology called cmPALS. This is a significant advance on existing PALS technology because it allows the modulator to make large movements. This means you can use shorter measurements and apply lower electric fields, reducing the effects of electrode fouling and deterioration.

Zeta potential specifications

Measuring range	-600 mV to +600 mV
Size range	3.8 nm to 100 μ m
Sensitivity	1 mg/mL (lysozyme)
Max. sample concentration	40 % w/v
Sample volume	350 μ L
Max. sample conductivity	200 mS/cm
Measurement angles	15°

Transmittance

Transmittance is measured by detecting the amount of light that passes through the sample. The Litesizer™ 500 continuously measures the transmittance for every sample. The value is reported in real time and is displayed during operation.

Your benefit:

You obtain instant insight into the suitability of the sample for light-scattering measurements. In addition, this measurement enables the Litesizer™ 500 to select the best parameters for your sample (focus position, measuring angle, measurement duration).

Transmittance specifications

Size range	No limit
Measuring time	10 s
Min. sample volume	20 μ L

Instruments for Particle and Materials Characterization



Anton Paar

The SurPASS 3 system enables the Zeta Potential analysis of solid surfaces, powders, fibers, foils and porous materials. Real-life samples can be measured directly without the need for model surfaces.

Applications/Materials

- ▶ Zeta Potential of Solid Surfaces
- ▶ Filters/Membranes
- ▶ Medical Devices
- ▶ Biomaterials
- ▶ Semiconductors
- ▶ Contact Lenses
- ▶ Core Samples

Measurements


- ▶ Streaming Potential
- ▶ Streaming Current
- ▶ Isoelectric Point
- ▶ pH Scan
- ▶ Adsorption and Desorption Processes

Unique Features/Benefits

- ▶ Fully Automated pH Scan, ranging from pH of 2–12
- ▶ Direct analysis of solid surfaces, no need for model surface or reference particles
- ▶ Measuring time as low as 20 seconds
- ▶ Multiple samples holders
- ▶ Software alerts for sample holder recognition and pH calibration



SurPASS 3



Titanium shows negative surface charge at physiological pH and an IEP at pH 4.2 due to its amphoteric behavior. Adsorption of bovine serum albumin (BSA) on the titanium surface shifts the IEP to pH 4.7.

**Contact us to schedule an on-site demo
and see how you can save on both!**

Instruments for Particle and Materials Characterization



The Litesizer™ 500 measures nano and microparticle size, zeta potential, molecular mass and transmittance using Light Scattering technology with ingeniously simple software.

Applications/Materials


- ▶ Particle Size and Size Distribution
- ▶ Particle Zeta Potential
- ▶ Protein Aggregation
- ▶ Particle Dispersion
- ▶ Formulation Development
- ▶ Polymers, Metal and Metal Oxide Particles

Measurements


- ▶ Dynamic Light Scattering (DLS)
- ▶ Electrophoretic Light Scattering (ELS)
- ▶ Static Light Scattering (SLS)
- ▶ Zeta Potential
- ▶ Molecular Weight
- ▶ Transmittance

Unique Features/Benefits

- ▶ Analysis at 15, 90 and 175 degrees
- ▶ Advanced algorithms for particle size analysis
- ▶ Patented cmPALS technology
- ▶ Novel Omega-shaped cuvette
- ▶ Transmittance analysis
- ▶ 6 minute start up time
- ▶ Resistant to dust, vibration and temperature fluctuations



Litesizer™ 500



Size (radius) distribution of a tridisperse polystyrene latex sample containing particles of 22, 100 and 700 nm.

**Contact us to schedule an on-site demo
and see how you can save on both!**